


JOHNS HOPKINS  
BIOMEDICAL ENGINEERING


# SMART OR DEVICE

AUTOMATING THE FUTURE OF OPERATING ROOM SAFETY


## BACKGROUND




Cauterizers cut tissue and are a heat source during surgery



When supplied oxygen is high, cautery increases fire risk



Over **600** OR fires occur every year across the United States



## WHEN SHOULD FIOR BE USED?

FiOR can be implemented in any operating room where anesthesia and cautery are used frequently. Device stays active throughout procedures

## WHAT DOES FIOR DO?

FiOR automates OR safety

High oxygen supply detected in presence of cauterizer

FiOR shuts off cauterizer and notifies surgeon

Anesthesiologist reduces O<sub>2</sub> to safe level

## TEAM MEMBERS

Emily Guan  
William Sun  
Ritvik Puly

Aydin Turkey  
Noah Lu  
Trisha Karani  
Selena Kim


## PROJECT ADVISORS

- Dr. Miranda: Dept of Biomedical Engineering, Johns Hopkins University
- Dr. Mahoney: Wilmer Eye Institute, Johns Hopkins Hospital
- Niam Mohseni: Dept of Biomedical Engineering, Johns Hopkins University


## PROTOTYPE AND VISION




Current Prototype



FiOR Module



Electrosurgical Unit



Connection to Cauterizer

Autonomous

Reliable

Minimal

Affordable